



VC SIMPLY

Multifunctional mouth protector with replaceable protective shield

Article no.	VC4001
Product name	Multifunctional mouth protector with 3 replaceable protective shields
Manufacturer	Cangnan Smile Transparent Plastic Mask Factory
Product description	<ul style="list-style-type: none"> • Chin-secured for superb hold • Ultra-light and robust workmanship • Extremely comfortable to wear thanks to elastic ear-loops • Anti-fog coating • Ideal for use with loupes • Visor can easily be cleaned with water or disinfectant • Suitable for medical applications
Characteristics	<ul style="list-style-type: none"> • Visor made of 300µm-thick polycarbonate film • High tensile strength • Chemically and thermally stable • Dimensions of protective shield: 142 x 62 mm
Material	Shield: Polycarbonat
Application	Mouth protector to minimize the risk of droplet infection via liquids and other particles.
How to wear	The protective shield has a protective film on both sides. Remove the films before use. Insert the protective shield into the notch in the frame. Position the frame on the chin using the fixation. Pull the ear loops around the ears.



VC SIMPLY

Multifunctional mouth protector with replaceable protective shield

Requirements for use	Before use, check the protective shield for any defects such as rips or holes. Keep the protective shield away from open flames and heat sources.
Color	White
Storage	Don't store in direct sunlight, store dry in original packaging.
Disposal	After use, the protective shield can be disposed with the household waste.